

KLM Orthotic Labs System_Rx Shells Sizing Chart - Please use the System_Rx Shells Order Form to place orders

Shoe Size	M1-2 W3-4	M2-3 W4-5	M3-5 W5-7	M5-6 W7-8	M6-7 W8-9	M7-8 W9-10	M8-9 W10-11	M9-10 W11-12	M10-11 W12-13	M11-12 W13-14	M12-13 W14-15
Length	118	123	128	138	148	158	168	178	188	198	208
	Heel 40 X019 118-40-3-14	Heel 42 X020 123-42-3-14	Heel 38 X126 128-38-3-17	Heel 47 X022 138-47-3-14	Heel 47 X014 148-47-3-14	Heel 52 X005 158-52-3-14	Heel 52 X008 168-52-3-14	Heel 52 X018 178-52-3-14	Heel 57 X013 188-57-3-14	Heel 62 X015 198-62-3-14	Heel 67 X131 208-67-3-17
		X129 123-42-3-17	Heel 42 X026 128-42-3-14	X123 138-47-3-17	X118 148-47-3-17	X101 158-52-3-17	X108 168-52-3-17	X117 178-52-3-17	X113 188-57-3-17	X115 198-62-3-17	X222 208-67-3-20
		Heel 45 X025 123-45-3-14	X127 128-42-3-17	Heel 52 X024 138-52-3-14	Heel 52 X009 148-52-3-14	X201 158-52-3-20	X208 168-52-3-20	X215 178-52-3-20	X214 188-57-3-20	X212 198-62-3-20	X319 208-67-3-23
		X130 123-45-3-17	Heel 45 X021 128-45-3-14	X120 138-52-3-17	X111 148-52-3-17	X313 158-52-3-23	X314 168-52-3-23	X315 178-52-3-23	X310 188-57-3-23	X309 198-62-3-23	
			X122 128-45-3-17	X216 148-52-3-20	X105 158-57-3-14	Heel 57 X002 168-57-3-14	Heel 57 X003 178-57-3-14	Heel 57 X006 188-62-3-14	Heel 62 X007 198-67-3-14	Heel 67 X016 208-72-3-17	Heel 72 X125 208-72-3-20
			Heel 47 X023 128-47-3-14	X017 148-57-3-14	X205 158-57-3-17	X202 168-57-3-17	X206 178-57-3-17	X207 188-62-3-17	X210 198-67-3-17	X211 208-72-3-20	X318 208-72-3-23
			X119 128-47-3-17	X116 148-57-3-17	X311 158-57-3-20	X301 168-57-3-20	X306 178-57-3-20	X304 188-62-3-20	X308 198-67-3-20		
				X217 148-57-3-20	X011 158-62-3-14	X104 168-62-3-14	X103 178-62-3-14	X010 188-67-3-14	X121 198-72-3-17	X128 208-77-3-17	
					X213 158-62-3-17	X204 168-62-3-17	X203 178-62-3-17	X209 188-67-3-17	X316 198-72-3-20	X221 208-77-3-20	
					X312 158-62-3-20	X303 168-62-3-20	X302 178-62-3-20	X305 188-67-3-20			
					Heel 62 X011 158-62-3-23	Heel 62 X001 168-62-3-23	Heel 62 X004 178-62-3-23	Heel 67 X012 178-67-3-14	Heel 72 X124 188-72-3-17	Heel 72 X218 188-72-3-20	
							X110 178-67-3-17	X210 178-67-3-20	X317 188-72-3-23		
							X307 178-67-3-23				

For more information
Contact us at:
Phone: 800-556-3668
Fax: 800-556-3338
Email: cservice@klmlabs.com

Sizing Procedure

- 1) Determine length by shoe size
- 2) Determine heel width
- 3) Determine medial arch height

* All measurements are in millimeters

Choose Shell Stiffness/Color

Standard Rx Shell (Semi-Rigid) Blue, Black, White

Rx Rigid-Flex (In Between Semi-Rigid and Semi_Flexible) Brown

Rx Semi-Flex (Semi-Flexible) Flesh Tone

Rx Flex (Flexible) Clear

* Use a separate form for each Color/Rigidity

